

Code BT 001, 002, 003

MOHReg 5/2005/1



Sterilization By EtO gas.

Main Features

- Blood tapered spike with sharp piercing tip, and wide internal tract.
- 20ml drip chamber including cylindrical blood filter conforming to DIN 58360
- Precise roller clamp.
- T shape injection site to assure safety to patient and staff.
- 125 cm. glass clear tubing.
- Luer adaptor, or Threaded luer lock adaptor to secure the hermetic locking and stability of luer / needle hub.
- Ultra thin-wall double beveled siliconized G19 x 1.5" needle on request.

Supplies on Request

- Blood transfusion sets with twin compartment drip chambers assuring easy priming, constant flow and avoiding air bubbles.
- Blood transfusion set

integrated with IV set.

Packing

- Economic HDPE pack assuring sterile internal tract.
- Standard blister pack (medical grade paper, PA / PE laminate) assuring sterility to both internal and external tract.