

**Surgical Washing Set  
( Urology Y set )**

**Code(s) :** 1/17

**Main Features :**

- \* Designed to precisely control the flow of solution from GLYCINE containers to the urethral catheter, resect, or cystoscope.
- \* Wide 4.8x6.8 tubing to assure high flow.
- \* Stout trochar assures good connection to the solution administration port.
- \* Connection to the catheter, resector, or cystoscope through soft, wide, flexible and translucent TPE tube.
- \* Step-wise multi-click clamp to regulate flow precisely.

**Supplies On-Request :**

Any special design can be supplied in short time.

**Packing :**

- Economic HDPE pack assuring sterile internal tract.
- OR - Standard blister pack ( medical grade paper, PP/PE laminate ) assuring sterility to both internal and external tracts.

**Sterilization :** By EtO gas.